

4. LABORATORY TEST METHODOLOGY

This methodology aims at reproducing the washing and flotation steps of the recycling process at a small scale to determine the suitability of inks applied on packaging for the HDPE or PP recycling stream. The methodology described below shall be followed precisely and any modifications or problems must be noted during the testing phase.

An Evaluation Report compiling all the results obtained shall be prepared by the laboratory to explicit the results. Any remarks during following the Procedure shall be also noted down. An Annex to this document is also available on RecyClass website to summarize all the collected data and do the assessment

See below in Figure 1 a diagram where the flow of the methodology is described.

Figure1: Methodology Diagram

